

## How to Use Your Genuair Pressair Inhaler

1. Remove the cap from the device by pressing the arrows shown on each side of the cap.



2. Make sure there is nothing in the mouthpiece.

3. Hold the device with the mouthpiece toward you and the green button facing straight up.



4. Press the green button all the way down.

5. **RELEASE** the green button.



6. Before breathing in, make sure the coloured control window has turned to green.



7. Blow out away from the inhaler.

8. Seal your lips around the mouthpiece and breathe in **STRONGLY** and **DEEPLY**.

9. You should hear a “CLICK” when you are breathing in.



10. Remove the device from your mouth and hold your breath for up to 10 seconds.

11. Make sure the coloured control window has changed to red. This means that you received the full dose.

12. If the control window stays green, repeat steps 8-10.

13. Replace the cap on the inhaler.



14. The green button will lock in place when the device is empty.

15. Clean your inhaler by wiping the mouthpiece with a dry cloth only.

**DO NOT wash your inhaler.**

16. Store the medication between 15-30 degrees Celsius.

17. The medication comes in a protective pouch and must be used within 90 days of opening the pouch.